

## Materials and Methods

### 1. Seamless cloning MasterMix constructs plasmid

#### 1.1 DNA manipulation.

PCR was performed with 2×Taq PCR MasterMix (Solarbio, Beijing, China). Plasmid DNA was isolated with the TIANGEN Mini Plasmid Kit (TIANGEN, Beijing, China). DNA fragments were purified from agarose gels by using the Universal DNA Purification Kit (TIANGEN, Beijing, China) or the Monarch DNA Gel Extraction Kit (NEB). DNA sequencing and primer synthesis were carried out by Sangon Biotech (Shanghai, China).

#### 1.2 Plasmid construction.

- a. The plasmids and primers used in this study are listed in the supplemental material.
- b. To construct pTarget-luxAB, the coding regions of luxAB, p15A, and AmpR are amplified from the genome of Fg-1, a plasmid from 2021OUC-China iGEM and Cloned UpB\_4A3m with three different primer pairs luxAB-F/luxAB-R, p15A-R/p15A-F, AmpR-R/AmpR-F, respectively. The J23119 promoter is added to PCR products of luxAB by primers luxAB-F and J23119R, resulting in J&luxAB. Oligonucleotides are designed to contribute flanking homologous regions to adjacent DNA fragments of 20-25bp in length, resulting in p15A-homo and AmpR-Homo.
- c. A new pEvolVR plasmid expressing enCas9-PolI3M-TBD with gRNA targeting luxAB is constructed. Three fragments pEA, pEB and pEC with flanking homologous regions to adjacent DNA fragments of 20-25bp are amplified from pEvolVR-enCas9-PolI3M-TBD with primer pairs pEA-R/pEA-F, pEB-R/pEB-F, pEC-R/pEC-F, respectively. Double-stranded gRNA is gained after DNA annealing, and it is modified by the promoter with the primer pair——. Similarly, then homologous regions are added to produce J&gRNA-Homo.
- d. Two plasmids were both constructed using the Seamless cloning Master Mix provided by Sangon Biotech (Shanghai, China). The plasmids pTarget-luxAB and pEvolVR were cotransformed into E. coli strain DH5 $\alpha$ , following the specification steps provided.

### 2. Structural forecasting

The nucleotide sequences of gene luxA and luxB of fg-1 were translated into amino acid sequences from the first promoter. The structures of both luxA and luxB protein were predicted by phyre2 (<http://www.sbg.bio.ic.ac.uk/phyre2/html/page.cgi?id=index>) through homology modeling method. Further analysis of protein structure including the alignment and sites labeling were accomplished by PyMol.

### 3. Fluctuation Test

#### 3.1. Optimize the scale of training

- Start LB cultures of each E. coli strain and grow cultures overnight at 37°C with orbital shaking at 170 rpm.
- Prepare selective media.
- Dilute the bacteria until the culture does not produce colonies on the resistance plate.
- Determine the cell density of the overnight culture as precisely as possible, and the total number of cells at this time is N<sub>0</sub>.
- Depending on the measured bacterial density, the overnight culture was diluted to 30 ml with a

complete medium with different glucose concentrations (0.1%, 0.05%, 0.01%, 0.005%) so that each 30 ml solution contained a total of 100,000 cells.

- Measure the cell density of the four conditions to determine dilution.
- Divided into two 96-well plates (48 holes per condition), each with 30ul.
- Incubate two plates (37°C) at the right temperature and do not shake.
- Obtain cell density for ten of the 48 holes under each condition to ensure that cell counts vary with glucose concentration.
- The remaining holes add 70ul distilled water to the 30ul culture and mix well before applying the plate to prevent the bottom residue. The coated plate remains dry.
- Calculate the number of colonies of zero-level events and find an appropriate proportion of zero-level events (unmutated plates/total plates between 10% and 80%)

### 3.2. Fluctuation test

- Overnight culture the bacteria liquid.
- Break up the bacteria liquid and measure cell density.
- Dilute the bacteria into a fully cultured solution at the selected glucose concentration.

### 3.3. Apply a separate strain

- Confirm the dilution multiple, then remove a certain volume of the bacteria fluid from the diluted bacteria to coat the common medium, and then count the colonies (which can be used as the initial number of cells in each hole).
- Add the diluted bacteria to 4 96-well plates.
- Two of the 96-well plates were cultured, and the other two were coated in the same selected medium, where only mutated bacteria could grow, and then count the colonies (the number of bacteria that have mutated- "a")
- After the culture of the two 96-well plates, ten of the holes were extracted to measure the cultured bacterial density, and the total number of cultured bacteria ( $N_f$ ) was calculated. Then the ten holes of bacteria liquid coated into a complete medium, calculate the number of colonies, at this time the number of colonies is the actual number of active bacteria (b), so as to obtain the survival rate of bacteria viability.
- Spread all remaining bacteria into the same selected medium and count (number of mutant cells- "c").

### 3.4. Data analysis:

Calculate mutation rates:

$$u = \frac{-\ln(p_0)}{\text{viability} \times N_f}$$

$$p_0 = \frac{N_f - c}{N_f - a}$$

$$\text{viability} = \frac{b - a}{N_f - a}$$